


ISO 15841 test

- 3-point bending
- knife edge supports
(radii 0.05 ... 0.13 mm)
- 10 mm distance between
outer supports




10 mm



Adenta

Thermadent 12



35.5° C +/- 0.5°


© Ibe/Segner

© Ibe/Segner



Adenta

Thermadent 14



35.5° C +/- 0.5°


© Ibe/Segner

© Ibe/Segner



Adenta

Thermadent 16



35.5° C +/- 0.5°


© Ibe/Segner

© Ibe/Segner



Adenta

Thermadent 16x22



35.5° C +/- 0.5°


© Ibe/Segner

© Ibe/Segner



Adenta

Adenta Flexadent NiTi 14



35.5° C +/- 0.5°

© Ibe/Segner

© Ibe/Segner

